



Xi'an Jiaotong-Liverpool University
西交利物浦大学

Standard Operating Procedure for: Flash 2000 Elemental Analyzer

**Department of Health and Environmental
Sciences**

By Yili Cheng

Contents

1. Purpose.....	3
2. Scope.....	3
3. Procedure	3
3.1 Instrument start-up	3
3.2 Determination of the blank value	6
3.3 Guide to run analyses.....	8
3.4 Shut off instrument	9
Document information:	9

1. Purpose

The purpose of this SOP is to regulate the operation, usage and maintenance of flash 2000 elemental analyzer.

2. Scope

This SOP applies to all personnel using and managing flash 2000 elemental analyzer.


3. Procedure

3.1 Instrument start-up

- 3.1.1 Open the valve of He and O₂ gas cylinder (in ES436A). If the total pressure is less than 3MPa, contact Yili.
- 3.1.2 Open the valve of He and O₂ at the back of bench (fig.1). Make sure the gas supply pressure is 0.3 MPa.



Fig.1 Opening the Gas valve

- 3.1.3 Switch on analyzer  (at the back of the instrument) and PC.

- 3.1.4 Click the icons (fig.2) show below in turn.




Fig.2 EA icon

- 3.1.5 Click on NCHS (fig.3).






Fig.3 Selecting NCHS mode



In the main menu select **View> View Maintenance** (fig.4).

The screenshot shows the 'View Maintenance' menu option highlighted in the 'View' menu. The 'Maintenance Program' dialog box is open, displaying a bar chart with two bars, 'Left 1' and 'Left 2', both showing a green bar with a yellow top section. The y-axis for 'Left 1' ranges from 0 to 150, and for 'Left 2' from 0 to 300. The dialog also includes fields for 'Number of runs until next maintenance', 'Number of runs since last maintenance', 'Number of runs to warning message', and 'Not used'. The 'OK' and 'Cancel' buttons are at the bottom right.

3.1.8 Click  to edit instrument status. Carrier flow: 140ml/min; oxygen flow: 250ml/min; reference flow: 100ml/min. Furnace and detector should be off currently. Click **send** to transfer the operating parameters to the instrument. (fig.5)



3.1.9 Click , select the option *special functions and leak test* (fig.6).

The left screenshot shows the 'Special Functions' tab in the software interface. The 'Leak Test...' button is highlighted with a red dashed box. Below it is the 'Auto-Zero Gas Channels' button. The 'Control' section shows 'Disable Sampling' with a radio button.

The right screenshot shows the 'Leak Test' dialog box. The 'Leak Test Status' section includes checkboxes for 'Carrier gas outlet close' and 'Reference gas outlet close', and input fields for 'Leak test time' (0 sec), 'Carrier flow' (140 ml/min), and 'Reference flow' (100 ml/min). The 'Start' button is highlighted with a red dashed box. Other buttons include 'Stop' and 'Done'.

Fig.6 Leak test

3.1.10 Click , set all parameters as below. Click send and transfer the setting to instrument.(fig.7)





Fig.7 Turning on the furnace, oven and gas supply

3.1.11 Click  to view current elemental status, wait till the flows reach the set value.

3.1.12 Click , select the option **detector**, in the section **TCD**, select **filament on** and click “send” (fig.8).





Fig.8 Turning on the detector

3.1.13 After about 50 minutes the instrument furnace reach the temperature settings and the LED Ready



on the synoptic panel lights up.



3.1.14 Click , select the option **Detector**. In the section **TCD**, click the button **Auto-adjust level at 1000 μ V** (fig.9).





Fig.9 Auto-adjust the detector level

3.2 Determination of the blank value

3.2.1 To fill sample table



Click , enter sample name, filename and weight. Select type and standard name. A typical sample table shows below (fig.10).

File Edit sample Verify chrom. file Balance Help						
	A	Sample name	Filename	Type	Standard name	Weight (mg)
1		blank1	blank1	Blank		
2		blank2	blank2	Blank		
3						
4		bypass1	bypass1	Bypass		
5		standardA	standard1	Std	Sulphanilamide	3.567
6		standardB	standard2	Std	Sulphanilamide	2.664
7		standardC	standard3	Std	Sulphanilamide	2.802
8		TF-7.5-2	TF-7.5-2	Unk		4.647
9		TF-7.5-3	TF-7.5-3	Unk		4.721
10		TF-7.5-4	TF-7.5-4	Unk		4.418
11		TF-7.5-5	TF-7.5-5	Unk		4.329

Fig.10 Sample table

The standard we usually use for soil is cystine.

Bypass: a standard substance is analyzed to condition the instrument.

Note:

- Filename can't be repeated, or your result will be covered.
- It is recommended to add a blank line after two blanks so that the sequence will stop after the blank samples and you can check the result.
- Analyze a standard (as a sample) after every 20 samples.

The sample table can be reedited during your analysis.

You can right click the blank before **sample name** to set this sample as actual sample being acquired (fig.11) when the system makes the wrong order.




Fig.11 Setting as actual sample being acquired

3.2.2 Click **recalculation > reset calibration factors** to reset the blank value and clear the previous calibrations (fig.12).




Fig.12 Reset calibration factors

3.2.3 Throw the first sample into the furnace and click  to start sequence (fig.13).





Fig.13 sequence start

3.2.4 Throw the second blank into the furnace before the first blank is finished.



- 3.2.5 Click  and then click ***select result > area*** in the following window to check the blank value (fig.14). **Acceptable limits of blank value: N<3000, C<10000, H<20000, S=0.** Contact LFST if the data don't falls into the limited range.




Fig.14 Checking the blank value

3.3 Guide to run analyses

- 3.3.1 Load bypass, standards and samples on the tray and load the auto sampler tray. Two trays could be used together. The second tray had to be loaded at the beginning if you need to use it or the order maybe wrong when it comes to the second tray.

- 3.3.2 Rotate the bottom tray clockwise (to the left) for one step.

Note: there should always be one sample in the cavity (the waiting area) in advance before it will be analyzed. So either throw it into the cavity or rotate the plate manually. (When the plate is loaded and sequence is in process, next sample will plunge into the waiting area automatically at the beginning of analyzing the current sample.)



- 3.3.3 Click  to start sequence.

You can tick the box “***shut-off detector temperatures and gas***” (fig.15) before click “***start now***” if the experiment will be finished late in the night or it is not convenient for you to operate the instrument.




Fig.15 Automatically shut-off detector temperatures and gas when sequence is finished

- 3.3.4 Click **View> View calibration curve** to check the calibration (fig.16).







Fig.16 View calibration curve

3.3.5 Click  to check the result. Click  to export result.

3.4 Shut off instrument

3.4.1 Click , tick the box furnace, oven, flow and filament to enable the off condition. Click **send** to send the command to the instrument.(fig.17)




Fig.17 Shut-off detector temperatures and gas

Caution: TCD must be shut off when cut off gas flow.

3.4.2 Wait till the temperature is below 200°C. Close the software. Turn off PC, analyzer and gases.

Document information:

Document No:	SOP-EA-01	Version No:	0.1 (draft)
Drafted/ modified by:	Yili Cheng	Effective date:	
Comments:			